

How to Use create_mat_file

The create_mat_file function generates a .mat file that is compatible with the MNG toolbox.

It allows you to save **channel data** and **comments** with precise timing information.

Syntax

create_mat_file(filename, 'channel', data, timestamps_or_fs_start, name, units, ...
'comments', comment_strings, comment_timestamps)

- filename — name of the file to save (string, **without .mat extension**).
 - 'channel' — keyword to start defining a data channel.
 - data — numeric array of the channel's data points.
 - timestamps_or_fs_start —
 - **Option 1:** An array of timestamps, same length as data,
 - **Option 2:** A 2-element vector: [sampling_frequency, start_time].
 - name — channel name (character array).
 - units — channel units (character array).
 - 'comments' — keyword to start defining comments.
 - comment_strings — cell array of comment texts.
 - comment_timestamps — array of timestamps for each comment.
-

Key Features

- Add **multiple channels** by repeating the 'channel' block.
 - Add **multiple comments** tied to specific timestamps.
 - Automatically saves sampling information, channel metadata, and comments into the .mat file.
-

Example

```
% Create example data
```

```
data1 = sin(0:0.01:10); % Channel 1: sine wave
```

```
data2 = cos(0:0.01:10); % Channel 2: cosine wave
```

```
sampling_frequency = 100; % Hz
start_time = 0; % Start at 0 seconds

comments = {'Start of recording', 'Midpoint', 'End of recording'};
comment_times = [0, 5, 10];

% Create the .mat file
create_mat_file('example_data', ...
    'channel', data1, [sampling_frequency, start_time], 'Sine Wave', 'V', ...
    'channel', data2, [sampling_frequency, start_time], 'Cosine Wave', 'V', ...
    'comments', comments, comment_times);
```

This will create a file called example_data.mat containing:

- Two channels ("Sine Wave" and "Cosine Wave") with their data, sampling frequency, and units.
- Three comments at 0 s, 5 s, and 10 s.

Notes

- **Timestamps vs Sampling Frequency:** If you already have precise timestamps for each data point, pass them instead of [fs, start].
- **Multiple Channels:** You can define as many 'channel' sections as you need.
- **Optional Comments:** 'comments' are optional. If omitted, only channel data is saved.
- **File Saving:** The .mat file will be saved in the current working directory unless you provide a path in filename.